



CONTACT US




Tell: 98 21-860 45 861
Fax: 98 21-860 45 864



Address: No.22, 4th Alley, Ghanbarzadeh St,
Beheshti Ave, Tehran - 1533814701 - Iran



www.kavoshcom.com



Why do people need to use ECG Patch?

Electrocardiogram enables your doctor to monitor your heart condition at rest or during physical activities. However, abnormal heart rhythms are often sporadic, therefore they need to be captured as they occur in order for physicians to diagnose and give treatment accordingly.

The main purpose of ECG Patch is to record your heart rate and rhythm during a long period of time in order to have better chances of detecting a symptom ("event").



ECG Patch

W e a r a b l e E C G

Everyday life under your
physician's watch




How to use

1. Install ECG Patch app on your mobile and connect it to your device via bluetooth.
2. Connect two disposable electrodes to the device and attach ECG Patch on your chest after selecting the placement site as specified by physician protocol.
3. Go about your day to day life as usual, while your ECG will be recorded up to 15 days thanks to the battery life time.
4. The application sends your ECG recordings to the cloud server to be stored.

Mobile Application: Just a few taps to start:

- Displaying a live, medically accurate ECG
- Sending EKG data to the cloud server instantly
- Detecting and annotating arrhythmia
- Storing data in mobile phone in offline mode
- Easy to use for elderly users



Green Heart Monitoring Cloud Service: You no longer feel away from your physician!

- Heart's Rhythm & activity records available & accessible at any moment
- Providing a complete medical grade report of your heart condition during a long period
- Detecting & marking heart irregularities as they occur
- Enabling your physician to monitor your heart function, live at any moment
- Unlimited data storage